

**Table S3.** Sensitivity and specificity of the assays used in this study. True positive (TP) is the number of positive results for one target obtained by the array on the positive control materials, false positive (FP) is the number of positive results for one target obtained by the array on the control materials which don't carry the target, false negative (FN) is the number of negative results obtained by the array on the control materials which contain the target, true negative (TN) is the number of negative results obtained by the array on control materials which don't carry the target. Sensitivity is the rate of true positive results on the total control materials tested which contain the target (TP/(TP+FN)), specificity is the rate of the true negative results on the total control materials tested which don't carry the target (TN/(TN+FP)).

Target	True positive	False positive	False negative	True negative	Sensitivity (True positive rate)	Specificity (True negative rate)	Note
<i>aac(6')-Ib-cr</i>	8	0	0	49	1	1	
<i>aac(3)-IIa</i>	6	1	0	50	1	0.98	May produce some cross reactions with <i>aac(3)-IId</i>
<i>aac(3)-IV</i>	2	0	0	55	1	1	
<i>ant(2'')-Ia</i>	4	0	0	53	1	1	
<i>ant(3'')-Ia</i>	36	2	3	17	0.92	0.89	May produce some cross-reactions with <i>ant(2'')-Ia</i> ; Decreased sensitivity due to lack of detection of some <i>aadA5</i> variants
<i>aph(3')-Ia</i>	13	1	1	43	0.93	0.98	
<i>bla<sub>KPC2/3</sub></i>	2	0	0	57	1	1	
<i>bla<sub>NDM1/2</sub></i>	1	0	0	56	1	1	
<i>bla<sub>PER1</sub></i>	1	0	0	57	1	1	
<i>bla<sub>CMY-1</sub></i>	1	0	0	57	1	1	
<i>bla<sub>CMY-2</sub></i>	7	0	0	6	1	1	These primers/probe were tested on a smaller strain panel
<i>bla<sub>CMY-3</sub></i>	8	0	0	54	1	1	These primers/probe are specific for CMY-2 + CMY-3. They should be used in combination with the primers/probe set CMY-2 to

							determine if CMY-2 or CMY-3 is present.
<i>bla</i> CTX-M9-group	5	0	0	53	1	1	
<i>bla</i> CTX-M2-group	3	0	0	55	1	1	
<i>bla</i> CTX-M1-group	25	0	0	36	1	1	
<i>bla</i> CTX-M8-group	3	0	0	57	1	1	
<i>dfrA1</i>	19	0	0	38	1	1	
<i>dfrA12</i>	5	0	0	52	1	1	
<i>dfrA17</i>	4	0	0	53	1	1	
<i>floR</i>	6	2	1	48	0.86	0.96	
<i>mcr-1</i>	6	0	0	51	1	1	
<i>mcr-2</i>	1	0	0	57	1	1	
<i>mcr-3</i>	3	0	0	56	1	1	
<i>mcr-4</i>	1	0	0	57	1	1	
<i>mcr-5</i>	1	0	0	56	1	1	
<i>mcr-9</i>	3	0	0	54	1	1	
<i>bla</i> OXA-48	2	0	0	56	1	1	
<i>qnrB1</i>	2	0	0	55	1	1	
<i>qnrS</i>	3	0	0	54	1	1	
<i>bla</i> SHV	9	0	0	50	1	1	
<i>strA</i>	31	0	0	26	1	1	
<i>strB</i>	31	0	0	26	1	1	
<i>sul1</i>	32	0	0	26	1	1	
<i>sul2</i>	34	0	0	23	1	1	
<i>sul3</i>	8	0	1	48	0.89	1	
<i>bla</i> TEM	43	0	0	24	1	1	
<i>tetA</i>	28	0	0	29	1	1	

<i>tetB</i>	12	0	0	45	1	1	
<b>IncA/C</b>	4	0	0	54	1	1	
<b>IncFIA</b>	11	0	0	46	1	1	
<b>IncFIB</b>	28	0	6	25	0.82	1	
<b>IncFII</b>	26	0	5	26	0.84	1	
<b>IncHI1A</b>	2	0	0	56	1	1	
<b>IncHI1B</b>	3	0	1	54	0.75	1	
<b>IncHI2</b>	10	0	0	47	1	1	
<b>IncI2</b>	3	0	0	54	1	1	
<b>IncI-I1</b>	15	0	1	41	0.94	1	
<b>IncK/B/O/Z</b>	9	0	0	48	1	1	
<b>IncL/M</b>	2	0	0	55	1	1	
<b>IncN-1</b>	1	0	0	57	1	1	
<b>IncN-2</b>	1	0	0	57	1	1	
<b>IncQ-1</b>	14	0	0	43	1	1	
<b>IncR</b>	4	0	0	53	1	1	
<b>IncX4</b>	9	0	0	49	1	1	Specific for IncX4_1, does not include IncX4_2
<b>ISSEn4</b>	7	1	0	49	1	0.98	
<b>ISKpn14</b>	1	0	0	57	1	1	Only exact ISKpn14 sequences are targeted
<b>ISKpn19</b>	4	0	0	54	1	1	
<b>ISKpn26</b>	23	1	1	33	0.96	0.97	
<b>IS26</b>	42	8	0	8	1	0.5	
<b>IS4321</b>	20	0	0	37	1	1	Targets both IS4321-L and IS4321-R
<b>ISKpn27</b>	1	0	0	57	1	1	
<b>IS903 + IS903B</b>	21	0	2	34	0.91	1	

<b>ISCR3/14/27</b>	2	0	0	55	<b>1</b>	<b>1</b>	
<b>ISCR1</b>	9	0	0	48	<b>1</b>	<b>1</b>	
<b>ISEc33</b>	1	0	0	57	<b>1</b>	<b>1</b>	Only exact ISEc33 sequences are targeted
<b>ISAb<sub>a</sub>24</b>	2	0	0	56	<b>1</b>	<b>1</b>	
<b>ISAb<sub>a</sub>125</b>	1	0	0	56	<b>1</b>	<b>1</b>	
<b>ISAp<sub>l</sub>1</b>	2	0	0	55	<b>1</b>	<b>1</b>	
<b>ISEcp1</b>	27	0	0	30	<b>1</b>	<b>1</b>	
<b>IS1294</b>	1	0	0	56	<b>1</b>	<b>1</b>	
<b>IS6100</b>	14	0	0	43	<b>1</b>	<b>1</b>	
<b>intI1</b>	41	2	2	14	<b>0.95</b>	<b>0.88</b>	
<b>intI2</b>	9	0	0	48	<b>1</b>	<b>1</b>	
<b>intI3</b>	1	0	0	57	<b>1</b>	<b>1</b>	
<b>Tn2</b>	32	0	1	24	<b>0.97</b>	<b>1</b>	